




# NutriCoat™ Attachment Solution

A xeno-free solution that facilitates the attachment and spreading of human mesenchymal stem cells cultured in serum-free medium

**REF** 05-760-1-15 (1.5ml/vial)

 15-25°C (Ambient temperature)

## Instructions for Use

### Product Description

A clear, colourless to yellowish, sterile solution, based on human Fibrinogen. NutriCoat™ Attachment Solution was qualified for the attachment of mesenchymal stem cells from various sources when culturing in MSC NutriStem® XF Medium.

### Precaution and Disclaimer

- Not for human use; For research or further manufacturing use only.
- Do not use solution that is cloudy or has deposits.
- Do not use beyond the expiration date indicated on the product label.
- In order to maintain the sterility of the product, it is necessary to use a strict aseptic technique and use sterile materials and tools only.

### Features

- Xeno-free
- Sterile by membrane filtration
- Easy-to-use
- Supplied at ambient temperature as a ready-to-use stock solution
- Design for use with serum-free, xeno-free culture system
- Enables cells dissociation using recombinant trypsin solution (e.g. 03-078-1)

### Instructions for Use

#### Pre-coating procedure

1. Dilute the solution 1:500 in saline and gently mix using a pipette.
2. Coat the culture dish with the 1:500 solution. The volume should be adequate for covering the surface of the culture ware. Use table 1 for the recommended volumes. (e.g. per 1 well of 12w/p add 2µl of NutriCoat™ in 1 ml saline)
3. Gently agitate and incubate for 1 hour in a humidified CO<sub>2</sub> incubator (37°C).

**Note:** longer incubation of coated culture dish (up to 24 hours) will not affect cell attachment and spreading.

4. Wash step is not required.

**Note:** one wash using Dulbecco's PBS (without Ca & Mg) is an option.

5. Prior to cells seeding, carefully discard the solution and add a culture medium.

**Note: It is critical that the coating does not dry out.**

**Table 1.** Recommended volume for the coating procedure

Culture ware	Surface area cm <sup>2</sup>	Volume of the 1:500 diluted solution
96-well	0.34	0.1 ml
24-well	1.9	0.5 ml
12-well	3.9	1 ml
6-well	9.6	2.5 ml
T25 Flask	25	6.5 ml
T75 Flask	75	19 ml

## Quality Control

Each lot is tested for sterility as well as its biological performance (the ability to promote the attachment and spreading of cells in serum-free medium).

## Quality Assurance

- For research use or for use as ancillary material in manufacturing cell, gene or tissue-based products.
- Composed of cGMP medical grade raw materials.
- Manufactured under ISO 13485 QMS and in compliance with applicable cGMP guidelines.
- Manufactured under controlled environments and processes in accordance with:
  1. ISO 13408 – Aseptic Processing of Health Care Products;
  2. ISO 14644 – Airborne Particulate Cleanliness Classes in Clean Rooms and Clean Zones.



Manufacturer:  
Biological Industries Israel Beit Haemek Ltd.  
Kibbutz Beit Haemek 2511500, Israel  
T.+972.4.9960595 F.+972.4.9968896  
Email: info@bioind.com

## Product Label Symbols



Indicates the manufacturer's catalogue number so that the product can be identified.



Indicates the manufacturer's batch code so that the batch or lot can be identified. Note: Synonyms for batch code are lot number and batch number.



Indicates the date after which the product is not to be used.



Indicates the temperature limits to which the product can be safely exposed.



Indicates a product that has been manufactured using accepted aseptic techniques.



Indicates the need for the user to consult the instructions for use.

## Related Products

Product	Cat. No.
MSC NutriStem© XF Basal Medium and Supplement Mix	05-200-1 05-201-1
MSC NutriStem© ACF Basal Medium and Supplement Mix	05-203-1 05-204-1
Recombinant Trypsin Solution	03-078-1
Recombinant Trypsin + EDTA Solution	03-079-1
NutriFreez™ D10 Cryopreservation Medium	05-713-1
Dulbecco's PBS (without Ca & Mg)	02-023-1
Saline - Sodium Chloride 0.9% Solution	03-043-1