# **RPMI 1640 product information**

PI-C3011 V1.0

## [Product Name]

Ival

Name: RPMI1640, with L-Glutamine, SLE Grade Cat. No: C3011-0500、C3011-0100 Specifications: 500ml、100ml

#### [Product Description]

Roswell Park Memorial Institute (RPMI) media are a series of media developed by Moore et al for the culture of human normal and neoplastic cells in vitro. RPMI 1640 is the most commonly used medium in the series. A modification of McCoy's 5A medium, the medium was specifically designed to support the growth of human lymphoblastoid cells in suspension culture. Presently the medium is extensively used for a wide range of anchorage dependant The medium needs to cell lines. be supplemented with 5-20% fetal bovine serum. The medium is also known to support growth of cells in the absence of serum.

It is super low endotoxin SLETM grade , Endotoxin  $\leq$  0.01EU/ml by LAL method. The production process and quality management system meet the requirements of cGMP and has passed ISO13485 certification.

Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

## [Composition]

# Ingredients

Mg/L

INORGANIC SALTS

Calcium nitrate tetrahydrate

100.000

Magnesium sulphate anhydrous	48.840
Potassium chloride	400.000
Sodium chloride	6000.000
Sodium phosphate dibasic	800.000
anhydrous	
AMINO ACIDS	
Glycine	10.000
L-Arginine hydrochloride	241.000
L-Asparagine	50.000
L-Aspartic acid	20.000
L-Cystine dihydrochloride	65.200
L-Glutamic acid	20.000
L-Glutamine	300.000
L-Histidine hydrochloride	20.960
monohydrate	
L-Hydroxyproline	20.000
L-Isoleucine	50.000
L-Leucine	50.000
L-Lysine hydrochloride	40.000
L-Methionine	15.000
L-Phenylalanine	15.000
L-Proline	20.000
L-Serine	30.000
L-Threonine	20.000
L-Tryptophan	5.000
L-Tyrosine disodium salt	28.830
L-Valine	20.000

🤽 400-820-3979 🕓 021-5878-5545 🔀 info@xpbiomed.com

1 / 2

# **RPMI 1640 product information**

PI-C3011 V1.0

Vitamins	
Choline chloride	

Choline chloride	3.000
D-Biotin	0.200
D-Ca-Pantothenate	0.250
Folic acid	1.000
Niacinamide	1.000
Pyridoxine hydrochloride	1.000
Riboflavin	0.200
Thiamine hydrochloride	1.000
Vitamin B12	0.005
i-Inositol	35.000
p-Amino benzoic acid (PABA)	1.000
OTHERS	
D-Glucose	2000.000
Glutathione reduced	1.000
Phenol red sodium salt	5.300
Sodium bicarbonate	2000.000

## [Storage and Stability]

RPMI 1640 , with L-Glutamine should be kept 2-8°C. The product is light -sensitive and therefore should not be left in the light. When stored in the dark under ideal conditions, the product is stable until the expiry date. As with any other liquid media formulations, deterioration of liquid media may be recognized by any of the following characteristics, among others including: (a). Color Change;

(b). Presence of clumping/flocculent debris/ granulation/ particulates\ precipitates or sediments ;

(c). Insolubility;

(d). And/or decrease in expected performance parameters.

Any material described above should not be used and therefore discarded.

## [Procedure]

1.Take a bottle from the defined storage conditions at 2-8°C and read the label.
2.Wipe the outside of the bottle with a disinfectant solution such as 70% ethanol.
3.Using aseptic/sterile technique under a laminar-flow culture hood, work according to established protocols.
4.Antibiotics may be added if desired.

## [Quality Contro]

RPMI 1640, with L-Glutamine is tested for sterility, pH, osmolality and endotoxin concentrations. In addition, each batch is tested for cell growth using Vero.

## [Precaution and Disclaimer]

For research use only, not for clinical diagnosis and treatment.