



Product Profile

Product Name:	Earle's Balanced Salt Solution (EBSS) without Phenol Red
Product Catalog Number	02-011-1
Concentration:	1X
Unit Size Availability:	(A)500ml:(B)100ml
Formulation:	Clear Solution
Optimal Storage Conditions:	Room Temperature (15-30°C)
Stability: (Under Specified Handling & Storage)	Please Refer To Product Label

Important Note! Please read the MSDS and Product Profile carefully in their entirety before using this material for possible safety precautions and potential hazards.

Product Description:

Balanced salt solutions are for all intents and purposes, inorganic salt solutions that have since been modified and enriched with a diverse complex of salt compounds, along with D-Glucose with or without Phenol Red that segue into a final salt solution based upon application and technique.

EBSS is one of several balanced salt solutions that form the basis of many complex media formulations. They are utilized to maintain cells for the short term in a viable condition rather than to promote their growth. In addition, these salt solutions may also be used for short incubations or for washing cells by centrifugation. During these intervals, the cells maintenance requirements are such that osmotic balance and physiological pH at the forefront.

These variegated inorganic salt solutions have been developed in order to fulfill the basic cell requirements for five basic and essential ions including: calcium, magnesium, phosphate, potassium and sodium and therefore contain various amounts of CaCl_2 , KCl , MgSO_4 , NaCl , NaHCO_3 , NaH_2PO_4 and other salts among others. The key constituents of salts are the ions which function in osmolality whereas others such as Calcium and Magnesium are known, among other functions, to serve as cofactors for and support cell attachment and aggregation. D-Glucose serves as a major carbon and energy source and Phenol Red may serve not only as a pH indicator but also may interfere, under some circumstances, with the growth of some cells at cloning densities in specific cell culture media. When added to this medium, Sodium Bicarbonate has an intimate relationship with CO_2 and plays a major role by helping to maintain optimal physiological pH.

The current role of a balanced salt solution in cell culture is multi-faceted and may be divided into four principal functions.

- Functions as a diluent, as an irrigating medium or transporting fluid while maintaining osmoregulation, the optimal and constant balance of osmotic pressure gradients between the intracellular and extracellular compartments.
- Provides cells with fluids and certain bulk inorganic ions essential for normal cell metabolism.
- When combined with a carbohydrate, such as $\text{C}_6\text{H}_{12}\text{O}_6$ (glucose), it provides a primary energy source for cell metabolism.
- When provided with a buffering system, it facilitates the maintenance of physiological pH within an acceptable range of 7.1-7.5.

Some Predominant Characteristics of EBSS without Sodium Bicarbonate and without Phenol Red include:

- ◆ Liquid Formulation
- ◆ With Sodium Bicarbonate (NaHCO_3)
- ◆ Without Phenol Red ($\text{C}_{19}\text{H}_{13}\text{NaO}_5\text{S}$) as pH indicator
- ◆ Sterile-Filtered (0.1 μ)
- ◆ Endotoxin-Tested

Instructions/Procedure:

The product should be stored at room temperature (15-30°C). The product should not be left in the light for prolonged periods as it is light-sensitive. When stored in the dark under ideal conditions, the product is stable until the expiry date.

- 1) Take a bottle from the proper storage conditions at room temperature (15-30°C) and read the label.
- 2) Ensure that the cap of the bottle is tight.
- 3) Gently swirl the solution in the bottle.
- 4) Wipe the outside of the bottle with a disinfectant solution such as 70% ethanol.
- 5) Using aseptic/sterile technique under a laminar-flow culture hood, work according to established protocols.

Quality Control:

Test	Specification
Appearance:	Clear Solution
Endotoxins:	Test and Record
Osmolality:	270-300 mOsm/Kg
pH :	7.1-7.5
Sterility:	Sterile

Auxiliary Products

Product Name	Catalog Number	Storage Temperature
Earle's Balanced Salt Solution without Phenol Red	02-010-1	Room Temperature (15-30°)
Earle's Balanced Salt Solution 10X Conc., without Sodium Bicarbonate	02-010-5	Room Temperature (15-30°)
Earle's Balance Salt Solution without Phenol Red, without Sodium Bicarbonate	02-011-5	Room Temperature (15-30°)
L-Glutamine Solution 29.2mg/ml in Saline	03-020-1	-20°C
L-Alanyl-L-Glutamine Solution(A Dipeptide Substitute)	03-022-1	-20°C
Sterile Culture-Grade Water	03-055-1	Room Temperature (15-30°)
Serum-Free Cell Freezing Medium	05-065-1	2-8°C
Note: For a list of antibiotics, serum and other reagents and supplements, please refer to our Product Catalog, Product Profiles, Product Guides and Internet Site.		

References:

- 1) Current Edition Merck Index
- 2) Biological Industries(BI) Specifications
- 3) Darling, D.C. and Morgan, S.J. Animal Cells: Culture and Media, New York: John Wiley & Sons,1994